Dental Group

3Shape Library

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Dental Group

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
<https://ipd2004.com/es/librerias-cad>

For internal use only. Do not distribute

Orders

- By date
 - Today
 - Last two days
 - Last week
 - Last two weeks
- By status
 - Created
 - Scanned (unprepared)
 - Scanned
 - Designed
 - Sent
 - Closed
- Other filters
 - Checked out
 - All open
 - My recent
- Search

Number	Creation date	External Lab	Delivery date	Customer	Items	Material	Status	Height	Expected delivery	Contact Person	Last modified	Manu...
2116146052_2...	14/09/2023 15:19:23	-	14/09/2023	1980192040	Crown 25, Scre...	Wax,Zirkon	Scanned (unprepared)	-/-/-	-		14/09/2023 15:27:24	1980...

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Home Page

System Settings

- System Settings
- Design options
- Services
- Auto workflow

DentalManager

- General
- Manufacturing Inbox
- ERP and processing time
- Order import
- Working days
- 3rd party applications
- Milling machines

Site Settings

- This site
- Dentists
- Labs
- Manufacturers
- Operators
- Countries
- Manufacturing processes

Tools

- Subscription Management (Dongle)
- Import/Export
- Download center
- Training Center
- 3Shape Communicate

Basic elements

- Materials
- Colors
- Margin line
- Attachments
- 2D Design overlays
- CAD blocks
- Press Multi sprues

Anatomy elements

- ScanIt library
- Anatomy and Pontic libraries
- Smile libraries
- Artificial teeth
- Crowns
- Crown Pontics
- Inlays
- Onlays/Veneers
- Temporary Crowns
- Temporary Pontics
- Temporaries on prepared model
- Tabletops

Frame elements

- Copings
- Frame Pontics
- Waxups
- Primary telescopes

Abutments

- Top cap libraries
- Implant systems

Post and Core

- Post and Core systems
- Post and Cores

Bridge elements

- Bars
- Connectors

Full dentures

- Design settings
- Gingivae
- Smile Libraries Chart

Removables

- Frames
- Connectors and Rests
- Wax profile strips
- Retentions
- Stippled waxes

Digital model

- Digital model design
- Articulator interfaces

Appliances

- Orthodontics Control Panel
- Positioning Guide

Seleccionar "Importar/Exportar"

Tools ▶ Import/Export

Tasks

- Export materials
- Import materials

Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".

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Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".

[Import materials](#)

History Of Import Materials

[Export materials](#)

[Clone system](#)

Open

Documents > 3Shape > Custom Interface

Name	Date modified	Type	Size
(OFF=00) 3Sh Zimmer TSV Custom Interface System Library 22-09-2022	1/16/2024 1:37 PM	DME File	16,834 KB
(OFF=25) 3Sh Straumann Bone Level Custom Interface System Library 21-09-2022	11/21/2022 1:41 PM	DME File	19,117 KB
(OFF=00) 3Sh Multi-unit Custom Interface System Library 04-02-2022	11/2/2022 2:59 PM	DME File	1,929 KB
(OFF=00) 3Sh Biomet 3i Externa Custom Interface System Library 04-02-2022	10/20/2022 12:59 PM	DME File	18,140 KB
(OFF=00) 3Sh Nobel Active Custom Interface System Library 04-02-2022	10/17/2022 2:49 PM	DME File	30,040 KB
(OFF=00) 3Sh Straumann Bone Level Custom Interface System Library 21-09-2022	9/26/2022 12:19 PM	DME File	19,116 KB

File name: (OFF=00) 3Sh Zimmer TSV Custom Interface System Library 22-09-2022

Dental System Materials Export

Open Cancel

Click "Open"

Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP -> Sites page and click "Download Materials".

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