



NEW  
COMPATIBILITIES  
AVAILABLE  
2024








Compatibility

# GALIMPLANT<sup>®</sup>

## 3A Serie MULTI-POSITION AESTHETIC<sup>®</sup> [ø3.8]





### INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Platform	Height (mm)		
ø3.8	10	IPD/3A-SR-00	IPD/3A-SR-01

Captive screw included. 1.20 Hex Screwdriver.




### Ti-BASE

Platform	Post height (mm)	Gingival height		
ø3.8	4.15	0.2	IPD/3A-IR-00/3D	IPD/3A-IR-01/3D
	6.5		IPD/3A-IR-02/3D	IPD/3A-IR-03/3D



### SCREW / TPA SCREW FOR ANGLED PROSTHESIS

Platform	Screwdriver	Seat	Reference	TiN Reference
ø3.8	Hex 1.30	Conical	IPD/3A-TR-00	IPD/3A-TR-00/TIN
			IPD/3A-TR-50	-
	TPA		IPD/3A-TR-51	-



### ANALOG / DIGITAL ANALOG

Platform	Reference	Digital reference
ø3.8	IPD/3A-AR-00	IPD/3A-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



## TEMPORARY ABUTMENT / OPEN TRAY COPING

Platform			Coping screw
ø3.8	IPD/3A-PR-06	IPD/3A-PR-07	IPD/3A-TR-00/C

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



## HEALING ABUTMENT

Platform	Height (mm)	Digital reference
ø3.8	4.5	IPD/3A-DR-04

IPD/3A-TR-00 screw not includes.





Compatibility



**CAMLOG<sup>®</sup>**

**JB Serie CONELOG<sup>®</sup>**

[ $\phi$ 3.3 /  $\phi$ 3.8 /  $\phi$ 4.3]





**INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY**

Platform	Height (mm)		
$\phi$ 3.3	10	IPD/JB-ST-00	IPD/JB-ST-01
	15	IPD/JB-ST-02	-
$\phi$ 3.8	10	IPD/JB-SN-00	IPD/JB-SN-01
	15	IPD/JB-SN-02	-
$\phi$ 4.3	10	IPD/JB-SR-00	IPD/JB-SR-01
	15	IPD/JB-SR-02	-



Captive screw included. 1.20 Hex Screwdriver.

**Ti-BASE  Custom Ti-Base System<sup>®</sup>**

Platform	Height (mm)		
$\phi$ 3.3	0.5	IPD/JB-IT-00/3D	IPD/JB-IT-01/3D
	1.5	IPD/JB-IT-02/3D	IPD/JB-IT-03/3D
$\phi$ 3.8	0.5	IPD/JB-IN-00/3D	IPD/JB-IN-01/3D
	1.5	IPD/JB-IN-02/3D	IPD/JB-IN-03/3D
$\phi$ 4.3	0.5	IPD/JB-IR-00/3D	IPD/JB-IR-01/3D
	1.5	IPD/JB-IR-02/3D	IPD/JB-IR-03/3D



**SCREW**

Platform	Screwdriver	Seat	Reference	TiN Reference
$\phi$ 3.3 $\phi$ 3.8 $\phi$ 4.3	Hex 1.25	Cónico	IPD/JB-TR-00 IPD/JB-TR-01	IPD/JB-TR-00/TiN

The reference IPD/JB-TR-01 is an abutment extractor screw.  
The IPD/JB-TR-01 Abutment Extractor Screw cannot be used if angled systems are used with the ConeLog<sup>®</sup> System.



**ANALOG / DIGITAL ANALOG**

Platform	Reference	Digital reference
$\phi$ 3.3	IPD/JB-AT-00	IPD/JB-AT-00/3D
$\phi$ 3.8	IPD/JB-AN-00	IPD/JB-AN-00/3D
$\phi$ 4.3	IPD/JB-AR-00	IPD/JB-AR-00/3D
Bottom Screw	-	IPD/KA-TA-00
Side Screw	-	IPD/KA-TA-01



## STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT



Straight, 17° & 30°

Platform	Height (mm)	Straight	17°	30°
ø3.8	1.5	IPD/JB-MN-01	-	-
	2.5	IPD/JB-MN-02	IPD/JB-MN-12	-
	3.5	-	-	IPD/JB-MN-33
ø4.3	1.5	IPD/JB-MR-01	-	-
	2.5	IPD/JB-MR-02	IPD/JB-MR-12	-
	3.5	-	-	IPD/JB-MR-33

Our straight Multi-units can be used for single and multiple restorations.



## TEMPORARY ABUTMENT / OPEN TRAY COPING

Platform			Coping screw
ø3.3	IPD/JB-PT-06	IPD/JB-PT-07	IPD/JB-TR-00/C
ø3.8	IPD/JB-PN-06	IPD/JB-PN-07	
ø4.3	IPD/JB-PR-06	IPD/JB-PR-07	

The same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



## PSD OVERDENTURE ABUTMENT

Platform	Height (mm)	Reference
ø3.8	2.0	IPD/JB-LN-02
	3.0	IPD/JB-LN-03
	4.0	IPD/JB-LN-04
ø4.3	2.0	IPD/JB-LR-02
	3.0	IPD/JB-LR-03
	4.0	IPD/JB-LR-04



## HEALING ABUTMENT

Platform	Height (mm)	Reference
ø3.3	2.0	IPD/JB-DT-02
	4.0	IPD/JB-DT-04
	6.0	IPD/JB-DT-06
ø3.8	2.0	IPD/JB-DN-02
	4.0	IPD/JB-DN-04
	6.0	IPD/JB-DN-06
ø4.3	2.0	IPD/JB-DR-02
	4.0	IPD/JB-DR-04
	6.0	IPD/JB-DR-06





C/ Rosa dels Vents, 9-15  
08338 Premià de Dalt (Barcelona)

Phone: (+34) 93 278 84 91  
info@ipd2004.com  
[www.ipd2004.com](http://www.ipd2004.com)

