

SCREW TECHNICAL SPECIFICATIONS

Description			Features							
REFERENCE	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE
IPD/AA-TN-00	Nobel Biocare® Branemark System®	Ø3.5	70°	2.25	M1.6x0.35	2.15	5.15	7.30	Unigrip	35
	Nobel Biocare® Nobel Active®	Ø3.7								
IPD/AA-TR-00	Nobel Biocare® Branemark System®	Ø4.1	75°	2.4	M2x0.4	1.5	5.85	7.35	Unigrip	35
	Nobel Biocare® Nobel Active®	Ø4.3								
IPD/AA-TR-01	Hexágono Externo (Plat. Regular)	Ø4.1	75°	2.5	M2x0.4	1.67	4.6	6.27	Hex. 1.20	20
IPD/AA-TR-02	Nobel Biocare® Branemark System®	Ø4.1	Plano	2.40	M2x0.4	1.45	5.9	7.35	Unigrip	35
	Nobel Biocare® Nobel Active®	Ø4.3								
IPD/AA-TW-00	Nobel Biocare® Branemark System®	Ø5.1	72°	3.40	M2.5x0.45	1.37	5.33	6.7	Unigrip	35
IPD/AB-TR-00	Multi-unit	Ø4.8	Plano	2.00	M1.4x0.3	1.35	2.1	3.45	Unigrip	15
IPD/AC-TN-00	Nobel Biocare® Nobel Replace® Select	Ø3.5	Plano	2.4	M1.8x0.35	1.8	6.5	8.3	Unigrip	35
IPD/AC-TR-00	Nobel Biocare® Nobel Replace® Select	Ø4.3 / Ø5.0	Plano	2.4	M2x0.4	1.55	8.35	9.9	Unigrip	35
IPD/AD-TT-00	Nobel Biocare® Nobel Active®	Ø3.0	70°	1.8	M1.4x0.3	2	5	7	Unigrip	15
IPD/BA-TR-00	Biomet® 3i® Externa™ 3i® (Plat. R)	Ø4.1	75°	2.4	M2x0.4	1.79	5.81	7.6	Hex. 1.20	20
IPD/BA-TR-01	Biomet® 3i® Externa™ 3i® (Plat. R)	Ø4.1	75°	2.4	M2x0.4	1.79	5.81	7.6	Hex. 1.25	35
IPD/BA-TR-02	Biomet® 3i® Externa™ 3i®	Ø3.25/ Ø4.1 / Ø5.0	Plano	2.4	M2x0.4	1.65	6.15	7.8	Cuad. 1.30	35
IPD/BB-TR-01	Biomet® 3i® Certain®	Ø3.4/ Ø4.1 / Ø5.0	Plano	2.00	M1.6x0.35	1.9	6.5	8.4	Hex. 1.20	20
IPD/BB-TR-02	Biomet® 3i® Certain®	Ø4.1 / Ø5.0	140°	2.30	M1.6x0.35	2.05	6.35	8.4	Hex. 1.20	20
IPD/CA-TR-00	Klockner® SK2-NK2	Ø4.2	90°	2.45	M1.8x0.35	2.39	3.56	5.95	Hex. 1.70	30
IPD/CA-TR-01	Klockner® SK2-NK2	Ø4.2	90°	2.45	M1.8x0.35	0.72	3.58	4.3	Cruz	30
IPD/CB-TR-00	Klockner® Essential Cone®	Ø4.5	90°	2.18	M1.8x0.35	2.3	4.9	7.2	Estrella	30
IPD/DA-TR-00	Straumann® Tissue Level (SynOcta® level)	Ø4.8	90°	2.4	M2x0.4	1.2	3.3	4.5	Torx	15
IPD/DA-TR-01	Straumann® Tissue Level (Implant level)	Ø4.6 / Ø6.0	30°	2.4	M2x0.4	3.03	3.67	6.7	Torx	35
IPD/DA-TR-02	Straumann® Tissue Level	Ø4.8	Plano	2.52	M2x0.4	1.45	5.9	7.35	Torx	35
IPD/DA-TR-03	Straumann® Tissue Level	Ø4.8	30°	2.57	M2x0.4	3.05	3.65	6.7	Torx	35
IPD/DB-TR-00	Straumann® Bone Level	Ø3.3 / Ø4.1	30°	2.1	M1.6x0.35	2.31	5.54	7.85	Torx	35
IPD/EA-TN-00	Astra® OsseoSpeed®	Ø3.0	120°	2.1	M1.4x0.3	2	5.5	7.5	Hex. 1.27	15
IPD/EA-TR-00	Astra OsseoSpeed®	Ø3.5 - Ø4.0	Plano	2.35	M1.6x0.35	1.95	6.3	8.25	Hex. 1.27	20
IPD/EA-TR-01	Astra® Cono 20°	Ø3.5-Ø4.0 / 4.5-Ø5.0	Plano	2.25	M1.4x0.3	1.85	2.95	4.8	Hex. 1.27	15
IPD/EA-TW-00	Astra OsseoSpeed®	Ø4.5 - Ø5.0	Plano	2.35	M2x0.4	1.9	8.4	10.3	Hex. 1.27	25
IPD/EB-TN-00	Astra® Evolution®	Ø3.6	60°	2.15	M1.6x0.35	2.12	5.33	7.45	Hex. 1.27	20
IPD/EB-TR-00	Astra® Evolution®	Ø4.2	100°	2.35	M1.8x0.35	1.92	5.70	7.62	Hex. 1.27	20
IPD/FA-TR-00	Zimmer® Screw Vent®	Ø3.5/ Ø4.5 / Ø5.7	140°	2.1	UNF 1-72	2.54	5.66	8.2	Hex. 1.25	30
IPD/FB-TR-00	Zimmer® SwissPlus®	Ø4.8	140°	2.3	M2x0.4	2.55	6.75	9.3	Hex. 1.25	30
IPD/FB-TR-01	Zimmer® SwissPlus®	Ø4.8	30°	2.4	M2x0.4	3.03	3.87	6.9	Hex. 1.25	30
IPD/GA-TN-00	Bti® Hexágono Externo	Ø3.5	Plano	2.4	M2x0.4	1.7	3.85	5.55	Hex. 1.20	20
IPD/GA-TW-00	Bti® Hexágono Externo	Ø5.5	Plano	2.55	M2x0.4	1.7	6.15	7.85	Hex. 1.20	20

Description			Features							
REFERENCE	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE
IPD/GB-TR-00	Bti® Interna Universal®	Universal® Ø4.1 / Ancha Ø5.5	Plano	2.4	M1.8x0.35	1.85	5.2	7.05	Hex. 1.20	20
IPD/HA-TN-00	Microdent® Microdent® System	Ø3.5	140°	2.5	M1.8x0.35	2.13	3.87	6	Hex. 1.20	20
IPD/HA-TR-00	Microdent® Microdent® System	Ø4.2 / Ø5.1	120°	2.7	M2x0.4	2.5	5.5	8	Hex. 1.75	30
IPD/IA-TR-00	Dentsply® Xive® Friadent®	Ø3.4 / Ø3.8 / Ø4.5	90°	2.05	M1.6x0.35	2.1	6.75	8.85	Hex. 1.20	25
IPD/IB-TR-00	Dentsply® Ankylos® Balance Base Narrow	C/X	60°	2.4	M1.6x0.35	1.66	3.04	4.7	Hex. 1.0	15
IPD/IB-TR-01	Dentsply® Ankylos®	C/X	60°	2.05	M1.8x0.35	2	6.15	8.15	Hex. 1.0	25
IPD/IB-TR-02	Dentsply® Ankylos®	C/X	120°	2.35	M1.8x0.35	2.65	8.45	11.1	Hex. 1.0	25
IPD/ID-TT-00	Osteoplus® Iridium™ Argentum®	Ø3.3 / Ø3.8 / Ø4.2 / Ø4.8	120°	2.40	M1.8x0.35	1.69	5.83	7.52	Hex. 1.25	25
IPD/JA-TR-00	Camlog® Camlog® Implant System	Ø3.3 / Ø3.8 / Ø4.3	90°	2.1	M1.6x0.35	2.35	8.1	10.45	Hex. 1.25	20
IPD/KA-CI-08	Hollow Screw	-	Plano	2.8	M2.65x0.3	0.7	1.4	2.1	Cruz IPD	10
IPD/LB-TR-00	BioHorizons® Tapered Internal	Ø3 / Ø3.5 / Ø4.5 / Ø5.7	90°	2.40	M1.8x0.35	2.42	4.83	7.25	Hex. 1.27	35
IPD/MB-TN-00	Sweden & Martina® Premium™ Kohno®	Ø3.3 / Ø3.8	120°	2.2	M1.8x0.35	1.78	6.32	8.1	Hex. 1.27	25
IPD/MB-TR-00	Sweden & Martina® Premium™ Kohno®	Ø4.25 / Ø5.0	120°	2.4	M2x0.4	1.78	6	7.78	Hex. 1.27	25
IPD/NA-TN-00	Phibo® TSH®	S2	90°	2.3	M1.8x0.35	2.1	5.55	7.65	Hex. 1.25	25
IPD/NA-TW-00	Phibo® TSH®	S5	120°	2.5	M2x0.4	2.5	5.5	8	Hex. 1.25	25
IPD/WB-TR-00	Megagen® Any Ridge®	Ø4.2	Plano	2.1	M1.8x0.35	2	6	8	Hex. 1.20	25
IPD/XA-TR-00	Dentium® Implantium®	Ø4.0 / Ø4.5 / Ø5.0	120°	2.3	M2x0.4	2	8.20	10.2	Hex. 1.27	20
IPD/XA-TR-02	Nobel Biocare® Branemark System®	Ø4.1	Plano	2.52	M2x0.4	1.45	5.9	7.35	Hex. 1.20	30
	Nobel Biocare® Nobel Active®	Ø4.3								
IPD/ZA-TR-00	Bredent® Blue Sky™	Ø3.5 / Ø4.5 / Ø5.5	120°	2.30	M1.8x0.35	2.2	8	10.2	Unigrip Bredent®	25

TPA SCREW TECHNICAL SPECIFICATIONS

REFERENCE	Description		Features							
	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE
IPD/AA-TN-50	Nobel Biocare® Branemark System®	Ø3.5	70°	2.3	M1.6x0.35	2.08	4.72	6.8	TPA	20
	Nobel Biocare® Nobel Active®	Ø3.5								
IPD/AA-TR-50	Nobel Biocare® Branemark System®	Ø4.1	75°	2.3	M2x0.4	2.02	3.78	5.8	TPA	20
	Sweden & Martina® Outlink®	Ø4.1								
	Phibo® TSH®	S3-S4								
IPD/AA-TR-52	Nobel Biocare® Branemark System®	Ø4.1	Plano	2.3	M2x0.4	1.4	4.4	5.8	TPA	25
	Biomet® 3i® Externa™ 3i®	Ø3.25 / Ø4.1 / Ø5.1								
	Bti® Hexágono externo	Ø4.1 / Ø5.5								
IPD/AA-TW-50	Nobel Biocare® Branemark System®	Ø5.1	72°	2.85	M2.5x0.45	1.61	4.8	6.41	TPA	25
IPD/AB-TR-50	Multi-unit	Ø4.8	Plano	2.3	M1.4x0.3	1.5	2.1	3.6	TPA	15
IPD/AC-TN-50	Nobel Biocare® Nobel Replace® Select	Ø3.5	Plano	2.25	M1.8x0.35	1.4	6.1	7.5	TPA	25
IPD/AC-TN-51	Nobel Biocare® Nobel Replace® Select	Ø3.5	Plano	2.25	M1.8x0.35	1.4	5.4	6.8	TPA	25
IPD/AC-TR-50	Nobel Biocare® Nobel Replace® Select	Ø4.3 / Ø5.0	Plano	2.3	M2x0.4	1.4	7.4	8.8	TPA	25
	Astra® OsseoSpeed®	Ø4.5 / Ø5.0								
IPD/AC-TR-51	Nobel Biocare® Nobel Replace® Select	Ø4.3 5 / Ø5.0	Plano	2.3	M2x0.4	1.3	6.6	7.9	TPA	25
IPD/AD-TR-50	Nobel Biocare® Nobel Active®	Ø4.3	75°	2.3	M2x0.4	2.02	4.98	7	TPA	25
IPD/AD-TT-50	Nobel Biocare® Nobel Active®	Ø3.0	70°	2.30	M1.4x0.3	1.99	5.01	7	TPA	15
IPD/BB-TR-51	Biomet® 3i® Certain®	Ø3.4/ Ø4.1 / Ø5.0	Plano	2.25	M1.6x0.35	1.35	6.25	7.6	TPA	20
IPD/CA-TR-50	Klockner® SK2-NK2	Ø4.2	90°	2.3	M1.8x0.35	1.50	3.00	4.5	TPA	25
IPD/CB-TR-50	Klockner® Essential Cone®	Ø4.5	90°	2.3	M1.8x0.35	1.65	4.85	6.5	TPA	25
IPD/DA-TR-50	Straumann® Tissue Level (SynOcta® level)	Ø4.8 / Ø6.0	90°	2.3	M2x0.4	1.65	3.32	4.97	TPA	25
IPD/DA-TR-51	Straumann® Tissue Level (Implant level)	Ø4.8 / Ø6.0	30°	2.3	M2x0.4	2.65	3.65	6.3	TPA	25
IPD/DB-TR-50	Straumann® Bone Level	Ø3.3 / Ø4.1	30°	2.3	M1.6x0.35	2.15	5.55	7.7	TPA	20
IPD/EA-TN-50	Astra® OsseoSpeed®	Ø3.0	120°	2.25	M1.4x0.35	1.65	5.45	7.1	TPA	15
IPD/EA-TR-50	Astra® OsseoSpeed®	Ø3.5 / Ø4.0	Plano	2.3	M1.6x0.35	1.50	6.30	7.8	TPA	20
IPD/EA-TR-51	Astra Uni 20° / 45°	All	Plano	2.3	M1.4x0.3	1.50	2.90	4.4	TPA	15
IPD/EB-TN-50	Astra® Evolution®	Ø3.6	60°	2.3	M1.6x0.35	2.19	5.11	7.3	TPA	20
IPD/EB-TR-50	Astra® Evolution®	Ø4.2	100°	2.3	M1.8x0.35	1.88	5.72	7.6	TPA	25
IPD/FA-TR-50	Zimmer® Screw Vent®	Ø3.5/ Ø4.5 / Ø5.7	140°	2.3	UNF 1-72	1.73	4.97	6.7	TPA	25
IPD/GA-TN-50	Bti® Externa Universal	Ø3.5	Plano	2.3	M2x0.4	1.5	3.6	5.1	TPA	25
IPD/GB-TR-50	Bti® Interna	Universal® Ø4.1 / Ancha Ø5.5	Plano	2.3	M1.8x0.35	1.5	5.0	6.5	TPA	25
IPD/HA-TN-50	Microdent® Microdent® System	Ø3.5	140°	2.3	M1.8x0.35	1.59	4.01	5.6	TPA	25
IPD/HA-TR-50	Microdent® Microdent® System	Ø4.2 / Ø5.1	120°	2.3	M2x0.4	1.85	4.15	6	TPA	25
IPD/IA-TR-50	Dentsply® Xive® Friadent®	Ø3.8	90°	2.3	M1.6x0.35	1.66	5.44	7.1	TPA	20
IPD/IB-TR-51	Dentsply® Ankylos®	C/X	60°	2.3	M1.8x0.35	1.82	6.88	8.7	TPA	25

 and IPD 2004 are registered trademarks of IMPLANT PROTESIS DENTAL 2004 SL. All other Trademarks mentioned are the property of their respective Owners and neither them nor the trademarks are associated to IMPLANT PROTESIS DENTAL 2004 SL.

<i>Description</i>			<i>Features</i>							
REFERENCE	COMPATIBILITY	PLATFORM	SEAT	HEAD DIAMETER	METRIC	HEAD (L1)	BODY (L2)	TOTAL (LT)	SCREW-DRIVER	TORQUE
IPD/LB-TR-50	BioHorizons® Tapered Internal	Ø3.0 / Ø3.5 / Ø4.5 / Ø5.7	90°	2.3	M1.8x0.35	1.89	4.61	6.5	TPA	25
IPD/MB-TN-50	Sweden & Martina® Premium™Kohno®	Ø3.3 / Ø3.8	120°	2.3	M1.8x0.35	1.52	5.08	6.6	TPA	25
IPD/MB-TR-50	Sweden & Martina® Premium™Kohno®	Ø4.25 / Ø5.0	120°	2.3	M2x0.4	1.55	4.85	6.4	TPA	25