



Compatibility

BTI[®]

GD Serie CORE*

[Ø3.5]

*Not compatible with 4.5mm implant length



SCANNING ABUTMENTS

INTRA-ORAL / DESKTOP / CONTACT PROBE SCAN BODY

Indication	Height (mm)	Ø 3.5
	10	IPD/GD-SN-00
	15	IPD/GD-SN-02
	10	IPD/GD-SN-01

Captive screw included. 1.20 Hex Screwdriver.



DEFINITIVE RESTORATIONS

Ti-BASE *ipd* Custom Ti-Base System[®]

Indication	Gingival Height	Ø 3.5
	0.5	IPD/GD-IN-00/3D
	1.5	IPD/GD-IN-02/3D
	2.5	IPD/GD-IN-04/3D
	0.5	IPD/GD-IN-01/3D
	1.5	IPD/GD-IN-03/3D
	2.5	IPD/GD-IN-05/3D



SCREWS

Surface	Screwdriver	Ø 3.5
Ti	1.20 Hex	IPD/GD-TN-00
TiN	1.20 Hex	IPD/GD-TN-00/TiN
TiN	TPA	IPD/GB-TR-50



STRAIGHT AND ANGLED MULTI-UNIT ABUTMENT

Angulation	Gingival Height	Ø 3.5
Straight	1.0	IPD/GD-MN-01
	2.0	IPD/GD-MN-02
	3.0	IPD/GD-MN-03
	4.0	IPD/GD-MN-04
17°	2.5	IPD/GD-MN-12
30°	3.5	IPD/GD-MN-33

Our straight Multi-units can be used for single and multiple restorations.



ANALOGS

Types	φ 3.5
Analog	IPD/GD-AN-00
Digital analog	IPD/GD-AN-00/3D
Bottom screw	IPD/KA-TA-00
Side screw	IPD/KA-TA-01



OVERDENTURE RESTORATIONS

PSD SYSTEM

Indication	Height (mm)	φ 3.5
	2.0	IPD/GD-LN-02
	3.0	IPD/GD-LN-03
	4.0	IPD/GD-LN-04



TEMPORARY RESTORATIONS

TEMPORARY ABUTMENT

Indication	φ 3.5
	IPD/GD-PN-06
	IPD/GD-PN-07
Long coping Screw	IPD/GD-TN-00/C

Same abutment can be used as a temporary abutment with the occlusal screw or as an impression coping using the coping screw.



HEALING ABUTMENT

HEALING ABUTMENT

Indication	Height (mm)	φ 3.5
	1.0	IPD/GD-DN-01
	2.0	IPD/GD-DN-02
	3.0	IPD/GD-DN-03



INSTRUMENTS AND ACCESSORIES

SCREWDRIVER

Indication	Lengths (mm)	φ 3.5
1.20 Hex	18mm	IPD/PH-20-18
	25mm	IPD/PH-20-25
	32mm	IPD/PH-20-32
TPA	18mm	IPD/KA-CT-18
	25mm	IPD/KA-CT-25
	32mm	IPD/KA-CT-32



ACCESSORIES

Indication	Lengths (mm)	φ 3.5
Clinical RA Shank Handle	16mm	IPD/MC-00-16
Laboratory RA Shank Handle	70mm	IPD/ML-00-70
RA Shank Ratchet	90mm	IPD/DI-00-90
Ratchet adapter	-	IPD/DI-KA-90
PSD and insert extractor tool	-	IPD/LL-PS-01

